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1	English introduction	4
1.1	General Information	4
1.2	Symbols used	5
1.3	Definitions (according to the norm EN ISO 17664)	6
2	Intended Use	7
2.1	Application area	7
3	Product description	7
3.1	Overview	7
3.2	Crank	8
3.3	Drill chuck	8
3.4	Luer-Lock-Connection	8
3.5	Fixing screw	8
3.6	Push button	8
4	Installation and use	9
4.1	First check	9
4.2	Warranty	9
4.3	Assembly/ disassembly	9
4.4	Operating instructions	9
5	Reprocessing	10
5.1	General requirements	10
5.2	Cleaning and disinfection	10
5.2.1	Pre-treatment	11
5.2.2	Machine cleaning and disinfection	11
5.3	Control	12
5.4	Maintenance	12
5.5	Packaging	13
5.6	Sterilization	13
5.7	Storage	13
5.8	Material durability	14

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1 English introduction


Thank you for purchasing this product.

You have acquired a product of German origin and in high quality.

1.1 General Information

Please read this manual carefully and keep it for future reference. Improper use may result in injury to the patient or cause damage to the product. When handling contaminated or biohazardous components / materials always observe general precautions.

Our delivery and mailing address is:

SMT Schilling Metalltechnik GmbH
Griesweg, 33
78570 Mühlheim
Germany

Phone: 07463 993090

Fax: 07463 9930912

E-Mail: <mailto:info@smt-metalltechnik.com>

Internet: www.smt-metalltechnik.de

Copyrights

All content, photos, texts and graphics are protected by copyright. They may not be copied, altered, reproduced or published without prior written permission


Reserve technical changes!

The illustrations and technical data used in this manual may differ slightly from the current state of the product.

1.2 Symbols used

Symbol	Signalwort	Gefahr	Folge
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Additional Symbols used:

Symbol	Beschreibung
	Symbol for products that are put into use according to the directive 93/42 EEC
	Cleaning and disinfection.
	NOTE: Additional hints or important information.
	Assembly.
	Disassembly.
	Manufacturer.
	Article/ Item number.
	Lot number.
NOT STERILE	Symbol for not sterile product.
	Follow instructions.

1.3 Definitions (according to the norm EN ISO 17664)

Chemicals for treatment

Recipes for solutions that are used for reprocessing

This includes for example detergents, tensides, cleaning solutions, disinfectants, enzyme cleaner and sterilizing media.

Cleaning

Removal of contamination from an object needed for further processing

Disinfection

Process used to reduce the number of viable microorganisms on a product to a previously specified level

Manual cleaning

Cleaning without using a cleaning/disinfection apparatus

Reprocessing

Activities, including cleaning, disinfection, or sterilization necessary to prepare a new or to reprocess a used medical device

Sterile

Free from any viable microorganisms

Sterilization

Validated procedure to kill off all forms of viable microorganisms on a device

Cleaning/Disinfection apparatus

An apparatus for the cleaning and disinfection of medical products and of other medical objects

Manufacturer

Organization or person responsible for the design, manufacture, packaging and labeling of a device before it is sold by its own name, irrespective of whether such operations are carried out by the person himself or by a third party authorized to do so.

Sterilizing device

A treatment chemical which has the ability to destroy microorganisms and viruses when used in the proper dilution / dosage and applied over the recommended exposure time.

Reprocessor

A Organization which is responsible for carrying out the actions that are necessary to prepare a new or used medical device for its intended use.

2 Intended Use

2.1 Application area

This product is used for clamping and hand-operating bone drills, thread cutters and similar products and may only be used by trained medical personnel, who can assess and control the risk to the patient through its use. Otherwise, the user assumes full responsibility.

DANGER!

The product will be delivered **non sterile**. Prior to the initial and each following use, the product has to be cleaned, disinfected, and sterilized as well as checked for visible irregularities and malfunctioning according to the indications given in this manual.

Guideline for medical devices (conformity declaration)
This product corresponds to the guideline 93/42/EEC.


3 Product description

This handheld drilling machine is used for clamping and hand-operating bone drills, thread cutters and similar products.

A huge advantage is that there are no oils or greases used in the complete device.

The transmission/ smooth-running gear is easy to clean and sterilize. The device has been designed so that the handheld drilling machine can be easily and quickly cleaned thanks to the intended cleaning connection. With about 547g total weight, the handheld drilling machine is very comfortable in the hand. The handle is anatomically adapted to the operator's hand and can be changed in seconds by a push-button from right to left-handed operation.

3.1 Overview

3.2 Crank

The crank operates the instrument. The drill chuck is driven by turning the crank. To disassemble the crank, it must be fixed and the fixing screw on the other side of the device must be removed anticlockwise. Now the crank can be pulled out of the drilling machine for cleaning or to change the position of the handle for left or right handed. The mounting is carried out in reverse order.

3.3 Drill chuck

The drill chuck is used to hold tools.

3.4 Luer-Lock-Connection


The Luer-Lock-Connection is used for cleaning/ drying the instrument.
See chapter 5.2 Cleaning and disinfection.

3.5 Fixing screw

It is needed to fix the crank and to position it on the desired side.

3.6 Push button

You have to press the push button to adjust the handle.

4 Installation and use

4.1 First check

Check the product immediately after receipt for eventual transport damages and completeness. Complaints can only be considered if the manufacturer or carrier is immediately notified. A damage report must then be submitted immediately to SMT Schilling Metalltechnik GmbH.

4.2 Warranty

Safety note: The responsibility for proper cleaning, disinfection and sterilization of instruments is incumbent upon the customer / product user. National arrangements, including limitations thereto, must be followed carefully.

SMT Schilling delivers exclusively tested and faultless products to their customers. All our products are designed and manufactured to meet the highest quality standards.

SMT Schilling as a manufacturer of products excludes from any warranty claims and assumes no liability for direct or consequential damages caused by:

- incorrect usage
- improper use, application or handling
- after usage on patients with Creutzfeldt-Jacob-Disease (CJD) or its variants, we refuse any responsibility for re-usage! We recommend destruction of instruments. A treatment and reuse also in accordance with the Directive RKI2 – done entirely on their own responsibility!
- improper treatment and sterilization
- improper maintenance and repairs
- failure to follow this instruction guide

Repairs may only be carried out by authorized **SMT Schilling** companies or individuals. Foreign-throughs result in the exclusion of all warranty claims.

4.3 Assembly/ disassembly

The drilling machine needs basically no assembly or disassembly.

However, when changing the position of the handle and for cleaning, the crank must be removed as described in chapter 5.2.

The plastic handle of the crank must be removed by a little effort, for optimal cleaning. After being postponed back in position it locks in again automatically.

WARNING!

The handle must always be positioned in the opposite direction of the crank, otherwise the drilling machine cannot be used and you can injure your hand or the device.

4.4 Operating instructions

1. When operating the crank, the instrument may only be held by the handle.
2. If changing the tool, the drill chuck may only be operated with the chuck key.
3. Always tighten the fixing screw completely.
4. For cleaning, the flushing hose must be securely attached to the Luer-Lock-Connection

DANGER!

The drilling machine consists of metal which can conduct electrical current and heat. Avoid contact with heat and voltage sources.

5 Reprocessing

5.1 General requirements

This instruction guide can't replace training, carefulness and user's knowledge of technology. We therefore put the relevant legislation, standards and recommendations (e.g. RKI2 or AKI1) assumed to be known (see "Standards / References") and therefore limit ourselves to be observed by the user for each instrument instructions and information which are of importance for our products. Reasons for these statements and the consequences for non-compliance risks are set out in the legislation and recommendations.

All products and components must be cleaned, disinfected and sterilized before each use; This also applies to the first-time use after delivery, since all products are delivered non-sterile (cleaning and disinfection after removing the packaging, sterilization after packaging). An effective cleaning and disinfection is an essential requirement for an effective sterilization.

After receipt, please check identity, completeness, integrity and function before bringing the instruments for treatment. Before every use of instruments, please prove them for cracks, deformations, damages and functional efficiency. Especially careful inspections are closures and all movable parts. Worn, corroded, deformed, porous or otherwise damaged instruments must be sorted out.

The user is responsible for the sterility of the product and that only sufficiently validated instrument- and product-specific procedures are used for cleaning, disinfection, or sterilization, that the apparatus used (disinfection apparatus, sterilization apparatus) are serviced and checked regularly, and that the validated parameters are respected in each cycle.

Please consider additionally the laws as well as the hygiene regulations of the medical practice or hospital. This applies specifically also to the differing requirements regarding effective prion inactivation.

5.2 Cleaning and disinfection

After use, the instrument can be cleaned with the following steps.

Connect the handheld drilling machine with a cleaning pistol (pressure max. 3 bar) over the existing luer connection. While flushing, turn the crank clockwise. As a result, the drilling machine is cleaned from the inside out, the rinsing liquid exits from the designated slots and at the shaft ends. In order to remove residual liquid from the device, compressed air is pressed into the interior via the luer connection by means of a compressed air pistol (oil-free, cleaned air), thus forcing the remaining liquid outwards. Again, the crank should be turned.

WARNING!

All products must be removed from the packaging before use. Other packaging, protective caps and protective foils etc. should be removed and disposed environmentally friendly.

NOTE

All products will be delivered non-sterile. Please clean, disinfect, and sterilize the products before the first use according to the instructions.

WARNING!

The crank, the crank handle and the fixing screw must be removed from the housing during mechanical cleaning/ disinfection (the clamped tool must be removed beforehand)

5.2.1 Pre-treatment

Directly after use (maximum 2 hours) major debris must be removed from the products.

This is done by the Luer-Lock-Connection for the interior. For the exterior use running water or a disinfection liquid; the disinfecting agent should be aldehyde-free (otherwise blood soiling will set), should have proven efficacy (e.g. DGHM- or FDA-licensed, or CE-marking), should be suitable to disinfect the product and compatible with the product (see section „material durability“).

ATTENTION!

Only use a soft brush or a clean soft cloth that are used only for this purpose for manual removal of debris, never use metal brushes or steel wool. Please note that the disinfectant used for the pre-treatment serves only to protect people and cannot replace the mandatory disinfection step that has to be performed subsequently after cleaning.

NOTE

Always observe the instruction for use from the cleaning and disinfection manufacturer.

WARNING!

Damage due to inappropriate cleaning! Use as much cleaning liquid in proper concentration as necessary and remove unnecessary cleaning liquid with a dry cloth.

WARNING!


Products made of aluminium alloys will be damaged by basic (pH > 9) cleaning liquids and solvents.

5.2.2 Machine cleaning and disinfection

Recommended procedure:

Equipment: Washer-disinfector (WD) from Miele with Vario TD program.
Neodisher MediClean forte from the company Dr. Weigert.

1. Place the instrument against movement in the WD, so that the spray directly hits the instrument. Clean parts only in disassembled and unstressed condition.
2. Start the Vario TD-program incl. thermal disinfection. The thermal disinfection has to occur in consideration of the A0-value and national provisions (prEN/ISO 15883).

Source and information to Miele: www.miele-professional.com

The process operates according to a fixed pattern:

- **Intensive pre-rinse** with cold water – 4 minutes

- **Emptying**

- **Cleaning** with neodisher MediClean forte at a temperature of 55°C – 10 minutes holding time.
Please note the corresponding dosing instructions of the manufacturer, in this case Dr. Weigert. (www.drweigert.com)

- **Emptying**

- **Neutralization** with warm water (> 40°C) – 3 minutes

- **Emptying**

- **Intermediate rinse** with warm water (> 40°C) – 2 minutes

- **Emptying**

- **Thermal disinfection** with deionized water (A0 value > 3000)

- **Drying** – 30 minutes. Remove the residual moisture from the instruments, especially in ports. To do this, use filtered compressed air, which does not lead to recontamination.

If there is still visible contamination (blood, tissue adhesions) after mechanical reprocessing, the cleaning and disinfection must be done again.

5.3 Control

- Instruments need to be cooled to room temperature!
- Assemble instruments for functional testing!
- Sort out damaged instruments and check success of cleaning disinfection (repeat if necessary), return with decontamination certificate to us.

5.4 Maintenance

Special maintenance is not required. The product should be sent back to the manufacturer or distributor if the action of moving parts or the clamping force is impaired, or in case of damage. Never use instrument oil!

NOTE

Repair must be carried out by SMT. If not respected, the warranty becomes void!

ATTENTION!

Please send back only reprocessed products to the distributor / manufacturer for repair (cleaned, disinfected, and sterilized).

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- Sterilization is done after a thorough cleaning. It does not replace the cleanliness of instruments to be sterilized.
 - Put the instruments into an approved packaging for steam sterilization (soft packaging, sterile Containers with appropriate storage).
 - Sterilization occurs according a validated steam sterilization process in fractioned Vacuum process at 132°C. Other methods are not allowed for this sterilization.
1. Fractionated pre-vacuum (3 times)
 2. Sterilization temperature 132°C plus tolerance according DIN EN ISO 17665-1 (or DIN EN 554)/ ANSI AAMI ISO 11134.
 3. Holding time in full cycle 5 minutes
 4. Drying time: 10 minutes

WARNING!

The crank, the crank handle and the fixing screw must be removed from the housing during disinfection (the clamped tool must be removed beforehand)

In order to avoid staining (silicate), or corrosion of the instrument, the steam must be free of ingredients. The recommended limit for particle contents in feed water and steam condensate is fixed by EN 285. Already infested instruments by corrosion must IMMEDIATELY be taken out of the Mineral processing cycle.

When sterilizing multiple instruments, the maximum capacity of the sterilizer must not be exceeded (see information by respective manufacturers).

Responsibility for cleaning, disinfection and sterilization of instruments rests solely with the operators of sterilization department or product users.

More detailed information on the reprocessing of medical devices can be found under www.a-k-i.org, the homepage of the "Working Group Instrument reprocessing". Here you can download for free the latest version of the RED brochure.

Never use the flash-sterilization procedure. Additionally do not use hot air sterilization, irradiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

5.7 Storage

- Dry, dust-proof, no external force, no major temperature fluctuations and not in the immediate vicinity of aggressive media.
- Useful in trays, containers, cabinets
- Otherwise, no special or additional requirements

5.8 Material durability

When choosing cleaning and disinfecting agents, please make sure that they do not contain the following:

- organic/ inorganic, or oxidizing acids (pH should not be less than 4, the use of neutral/enzymatic agents is recommended)
- bases (pH should not be higher than 9, the use of neutral/enzymatic agents is recommended)
- Solvents (alcohol, benzene, acetone,...)
- Phenol
- Chlorine, bromine, iodine
- Chlorine salts (in particular ammonium chloride compounds), chlorinated/halogenated hydrocarbons
- Oxidizing agents, peroxide, hypochlorite

Do not clean any of the products with metal brushes or steel wool.
The products must not be exposed at temperatures higher than 137°C (279 °F)!

5.9 Re-usability

The products can be re-used if they are not damaged and has been cleaned and sterilized according to the instructions provided. The use of a damaged or dirty product is in the responsibility of the user.

If a damaged product is used or if a product is re-used without having been cleaned and sterilized, then the company is not liable for any damage causes.
In case of disrespect all liability is excluded.

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